


Quinton Q3000 and Q4000 Tango Interface Notes

You can setup the Tango with your stress system using the following directions. You will need the ECG trigger cable specified below.

The ECG trigger cable is required for the Tango monitor to be able to take a BP measurement.



1. Plugging in the ECG trigger cable

Connect: the ECG trigger cable, SunTech part #91-0018-00 (1/8" phono jack to BNC)

From: the BNC External ECG connection on the rear panel of the Tango (BNC)

To: the QRS SYNC connection on the back of the stress system (1/8" phono jack)

NOTE: If there is another instrument already plugged into QRS SYNC (i.e. echo or pulse oximeter) then you can either unplug the existing instrument or use a Dual Mini Phone Jack Adapter to allow the Tango and the other instrument to use the QRS SYNC at the same time.

2. Setting up the Tango monitor

- Turn on the Tango using the switch at the rear top left corner of the Tango.
- Once the operating screen is displayed, press the *MENU/SELECT* button. Use the UP or DOWN buttons to select **Utilities** and press the menu button.
- Use the arrow buttons to select **Device** and press the *MENU/SELECT* button to edit. Choose **QUINTON** using the arrow buttons and press the select button to confirm the choice.
- Select **ECG trigger** and press the select button to edit. Choose **DIGITAL**↑ and press the select button to confirm the choice.
- Select **Exit** twice to return to the operating screen.