


Schiller CS 200 Cycle Interface

You can setup the Cycle with your stress system using the following directions. You will need the RS-232 cable specified below. It establishes communication between the monitor and your stress system. This cable allows your stress system to prompt the monitor when it needs a BP measurement. In addition, it allows the BP measurements taken by the Cycle to be transferred to your stress system display and reports.



RS-232 cable
Part# 91-0013-00



1. Connecting the RS-232 cable:

Connect: the RS-232 cable, SunTech part #91-0035-00 (9 pin female to 9 pin male).
From: the RS-232 connection on the connector panel of the Cycle (9 pin male).
To: BP port on the back of the stress system (9 pin female).

2. Setting up the Cycle monitor

- From the measurement view, press and hold the **SELECT** and **ADVANCE** buttons simultaneously for 3 seconds. The System Menu will be displayed.
- Press the **MENU** button once to proceed to the Communication Screen. Press **ADVANCE** until CS-200 is displayed.
- Press **MENU** twice to return to the measurement view.

4. Setting up the stress system:

- Go to **System > Settings > System Configuration > Other System**.
- Set the following parameters:
 - **QRS: On**
 - **QRS Trigger polarity: Pos**
 - **QRS Pulse length: 45 ms**
- Go to **System > Settings > System Configuration > Exercise**.
- Set the following parameter: **Bloodpressure unit: BP-200 or EBM-503**.