


# Cambridge Heart CH 2000 Cycle Interface

You can setup the Cycle with your stress system using the following directions. You will need the RS-232 cable specified below. It establishes communication between the monitor and your stress system. This cable allows your stress system to prompt the monitor when it needs a BP measurement. In addition, it allows the BP measurements taken by the Cycle to be transferred to your stress system display and reports.



RS-232 cable  
Part# 91-0065-00



## 1. Plugging in the interface cable

Connect: the cable, SunTech part #91-0065-00 (9 pin female & BNC to 9 pin female)

From: the RS-232 connection on the connector panel of the Cycle (9 pin male)

Note: the adjacent BNC connector will be left unconnected

To: the NIBP port on the back of the CH2000 (9 pin female)

## 2. Setting up the Cycle monitor

a. From the measurement view, press and hold the SELECT and ADVANCE buttons simultaneously for 3 seconds. The System Menu will be displayed.

b. Press the MENU button once to proceed to the Communication Screen. Press ADVANCE until CH2000 is displayed.

c. Press MENU twice to return to the measurement view.

## 3. Setting up the stress system (Note: requires software version 2.0.3 or greater)

a. From the EXERCISE SYSTEM CH 2000 screen, select **SYSTEM SETUP**.

b. Select **DEVICES I/O**.

c. Select **NIBP**.

d. Select **SunTech Tango**.

e. Select **OK**. If a window opens and asks if you would like to back up system, select **Don't Backup**.

- BP Prompts: Manual blood pressure readings can be initiated on the keyboard by pressing BP/Event on the numeric keypad.

- For assistance with the protocols, call Cambridge Heart Service at 800-CAM-WAVE.