


Schiller AT 10 and AT 60 Cycle Interface

You can setup the Cycle with your stress system using the following directions. You will need the RS-232 cable specified below. It establishes communication between the monitor and your stress system. This cable allows your stress system to prompt the monitor when it needs a BP measurement. In addition, it allows the BP measurements taken by the Cycle to be transferred to your stress system display and reports.



RS-232 cable
Part# 91-0035-00



1. Connecting the RS-232 cable:

Connect: the RS-232 cable, SunTech part #91-0035-00 (9 pin female to 9 pin male).
From: the RS-232 connection on the connector panel of the Cycle (9 pin male).
To: BP port on the back of the stress system (9 pin female).

2. Setting up the Cycle monitor

- From the measurement view, press and hold the **SELECT** and **ADVANCE** buttons simultaneously for 3 seconds. The System Menu will be displayed.
- Press the **MENU** button once to proceed to the Communication Screen. Press **ADVANCE** until *Quest* is displayed.
- Press **MENU** twice to return to the measurement view.

3. Setting up the stress system:

- From the main monitor screen, select **NIBP**. This will display the BP menu.
- Select **3 NIBP Setup** by pressing “3” on the keyboard. This will display the currently selected BP monitor, which should read **1 NIBP: EBM 502/3**. If this is not displayed press **ENTER** and select this device from the options.
- The “**FNCT**” key will return you to the main monitor screen.