



# SunTech Oscar 2

## Technical Specifications

### The Oscar 2 System Includes...

Oscar 2 monitor  
 AccuWin Pro™ 4 (AWP4) software (Download)  
 2 Orbit™ ABPM cuffs (Adult and Adult Plus) Carrying pouch & patient belt  
 USB interface cable  
 User guide  
 Standard 2 year warranty

### Oscar 2 System Monitors

99-0133-00 Oscar 2 Monitor - Standard  
 99-0133-01 Oscar 2 Monitor - Bluetooth  
 99-0133-02 Oscar 2 Monitor - w/ SphygmoCor Inside  
 99-0133-03 Oscar 2 Monitor - w/ SphygmoCor Inside - Bluetooth  
 99-0133-10 Oscar 2 Monitor - Standard - Without AWP4  
 99-0133-11 Oscar 2 Monitor - Bluetooth - Without AWP4  
 99-0133-12 Oscar 2 Monitor - w/ SphygmoCor Inside - Without AWP4  
 99-0133-13 Oscar 2 Monitor - w/ SphygmoCor Inside - BT - Without AWP4  
 98-0056-00 Oscar 2™ PowerPack Accessory Kit

### Orbit™ Comfort Blood Pressure Cuffs

98-0239-01 #1 Orbit™ ABPM Cuff (18-26 cm)  
 98-0239-02 #2 Orbit™ ABPM Cuff (26-34 cm)  
 98-0239-03 #3 Orbit™ ABPM Cuff (32-44 cm)  
 98-0239-04 #4 Orbit™ ABPM Cuff (42-55 cm)

### Straps and Belts

98-0032-07 Oscar 2™ Pouch  
 98-0037-00 Oscar 2™ Belt

### Extended Warranty Program

99-0027-00 One Year Extended Warranty  
 99-0027-03 Three Years Extended Warranty

### Extended Warranty Terms and Conditions

The Extended Warranty must be purchased either:

- At the time of the original sale
- Prior to the end of the initial warranty period
- Prior to the end of the purchased extended warranty period





## True24 App for Bluetooth Systems

Available in the Apple App Store



### General Specifications

**Validations:** BHS (A/A), ESH, ANSI/AAMI SP10:2002

**BP Technique:** Oscillometry with step deflation

**BP Range:** 25-260 mmHg

**HR Range:** 40-200 bpm

**Sampling Periods:** 3 adjustable (awake, asleep, and special)

**Time Intervals:** 5, 10, 15, 20, 30, 45, 60, 90 and 120 min.

**Event Key:** Program to allow additional patient readings

**Data Storage:** Flash memory stores up to 250 readings

**Power Source:** Two AA batteries

**Weight:** 284 grams, including batteries

**Dimensions:** 12 x 7 x 3 cm

**PC Requirements:** USB port, Windows® XP or later

**App Requirements:** Apple iPhone

### CPT Codes for ABPM Reimbursement

**CPT code: 93784:** Ambulatory blood pressure monitoring utilizing a system such as magnetic tape and/or computer disk, for 24 hours or longer; including recording, scanning analysis, interpretation, and report.

**ICD-10 code: R03.0:** White-coat hypertension. Elevated blood pressure reading without diagnosis of hypertension.

*\*Private insurance carriers may reimburse for additional indications under relevant ICD-10 codes.*

507 Airport Boulevard, Suite 117  
Morrisville, NC 27560-8200 USA

Tel: + 1.919.654.2300  
1.800.421.8626

Fax: + 1.919.654.2301  
SunTechMed.com



EMERGO Europe  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands

105 HuanGuan South Road, Suite 15 2~3/F

DaHe Community Guanlan,  
LongHua District, Shenzhen  
GuangDong PRC 518110

Tel.: + 86.755.29588810  
+ 86.755.29588986 (Sales)  
+ 86.755.29588665 (Service)

Fax: + 86.755.29588829  
SunTechMed.com.CN