



OPTIONAL REPLACEMENT CONNECTORS			
<i>Photo</i>	<i>Description</i>	<i>SunTech Part Number</i>	<i>Monitors</i>
Screw Type (Threaded)			
	Plastic Male Screw-Type (Threaded)	64-0076-01	SunTech 247 GE Critikon Dinamap Criticare 507 Welch Allyn Altas Welch Allyn Vital Signs Welch Allyn Spot Spacelabs Nihon Kohden (Lifescop L)
Slip Luer			
	Metal Male Luer Slip	64-0047-05	
	Plastic Male Luer Slip	64-0123-00	CASMED 740/750
	Plastic Female Luer Slip	64-0125-00	Baum
Quick Connect			
	Plastic Male Quick Connect	64-0047-03	Nihon Kohden Datex-Ohmeda*
	Plastic Female Quick Connect	64-0047-02	SunTech Tango+ SunTech Oscar 2 GE/Marquette two-tube Datex-Ohmeda* Protocol Propaq Welch Allyn Propaq
Bayonet (Rectus)			
	Metal Bayonet (Rectus)	64-0047-04	Philips/HP Siemens Colin Spacelabs Datascope (Accurtor Plus & Passport 2)
	Plastic Bayonet (Rectus)	64-0122-00	Philips/HP Siemens Colin

	Plastic Bayonet (Rectus)	64-0122-00	Spacelabs Datascope (Accurtor Plus & Passport 2)
Inline			
	Plastic Female Inline	64-0124-00	Marquette

*Datex-Ohmeda uses one male and one female quick connect on its two-tube monitors