

Clinical Grade  
Blood Pressure Cuffs

**SunTech**® *All Purpose Cuffs*  
CLINICAL GRADE BP CUFF

Cuffs Designed To Meet All Your Expectations



+ Durable enough for multiple measurements on multiple patients

+ Easy to clean and maintain

+ Suitable for the widest variety of clinical applications

+ Designed with patient comfort in mind

SunTech Medical®, a company with a twenty-five year history of research in blood pressure technology brings you the SunTech® All Purpose BP Cuff. Ranging from Infant to Thigh, including long sizes, SunTech's color-coded, All Purpose BP Cuffs are designed with a tapered end and an easy to clean high quality nylon that provides an outstanding combination of durability, patient comfort and longevity. SunTech's All Purpose BP Cuffs are durable enough to sustain multiple blood pressure measurements on multiple patients in the widest variety of clinical applications and environments.




Two-Year Warranty on the SunTech® All Purpose BP Cuff

**+ Features & Benefits + Description**

- Durability:** Long-lasting, high quality nylon material, with a tapered end, provides a durable yet comfortable multi-patient cuff solution.
- Cleaning/Disinfecting:** The same high quality nylon material combined with a removable bladder, makes the SunTech All Purpose BP Cuffs easy to clean, disinfect or machine wash and provides a faster drying time.
- Ease of Use:** Color-coding and interior range markings ensure proper cuff sizing and fit.
- Reliable:** 100% leak tested for trust, value, and reliability.
- Regulatory Requirements:** Designed to exceed AHA, AAMI-SP10, and MDD guidelines for results you can trust.
- Flexibility:** Single or double tube design with multiple connectors provides needed flexibility.
- Latex-Free/PVC-Free:** Designed to be 100% latex-free and PVC-free to maximize patient safety and minimize adverse reactions.
- Customizable:** SunTech provides private labeling solutions to meet your individual branding and marketing needs

**+ Ranges**



**+ Symbols**

- Arm circumference.
- Product in compliance with Council Directive 93/42/EEC Medical Devices Directive.
- Product contains no latex material.
- Product contains no PVC (polyvinyl chloride) material.
- Artery symbol and arrow should be placed over brachial or femoral artery.
- Index Line.
- Cuff index line must fall within range markings for an accurate measurement.



SunTech Medical, Inc.  
507 Airport Boulevard, Suite 117  
Morrisville, NC 27560-8200 USA  
Tel: + 1.919.654.2300  
1.800.421.8626  
Fax: + 1.919.654.2301

SunTech Medical, Ltd.  
Oakfield Estate  
Eynsham, Oxfordshire  
England  
Tel: +44 (0) 1865.884.234  
Fax: +44 (0) 1865.884.235

SunTech Medical, Ltd.  
Level 19, Two, International Finance Centre  
8, Finance Street, Central  
Hong Kong  
Tel: +852.2251.1949  
Fax: +852.2251.1950

