

One-Tube Design Specifications

Cuff Size (Print Color)	Range (cm)	Bladder Length x Width (mm)	Overall Length x Width (mm)	Female Quick Connect	Male Quick Connect	Threaded (Screw)	Female Slip Luer (Plastic)	Male Bayonet (Metal)	Male Bayonet (Plastic)	Male Slip Luer (Metal)	Male Slip Luer (Plastic)	No Connector
Infant (Orange)	8-13	110 x 38	230 x 63	98-0080-01	98-0084-31	98-0084-21	98-0080-31	98-0084-41	98-0084-71	98-0084-61	98-0084-81	98-0084-11
Child (Green)	12-19	150 x 68	315 x 93	98-0080-02	98-0084-32	98-0084-22	98-0080-32	98-0084-42	98-0084-72	98-0084-62	98-0084-82	98-0084-12
Child Long (Green)	12-19	150 x 68	445 x 93	98-0080-03	98-0084-33	98-0084-23	98-0080-33	98-0084-43	98-0084-73	98-0084-63	98-0084-83	98-0084-13
Small Adult (Royal Blue)	17-25	193 x 93	440 x 118	98-0080-04	98-0084-34	98-0084-24	98-0080-34	98-0084-44	98-0084-74	98-0084-64	98-0084-84	98-0084-14
Small Adult Long (Royal Blue)	17-25	193 x 93	515 x 118	98-0080-05	98-0084-35	98-0084-25	98-0080-35	98-0084-45	98-0084-75	98-0084-65	98-0084-85	98-0084-15
Adult (Navy Blue)	23-33	265 x 125	505 x 150	98-0080-06	98-0084-36	98-0084-26	98-0080-36	98-0084-46	98-0084-76	98-0084-66	98-0084-86	98-0084-16
Adult Long (Navy Blue)	23-33	265 x 125	645 x 150	98-0080-07	98-0084-37	98-0084-27	98-0080-37	98-0084-47	98-0084-77	98-0084-67	98-0084-87	98-0084-17
Large Adult (Burgundy)	31-40	320 x 158	645 x 183	98-0080-08	98-0084-38	98-0084-28	98-0080-38	98-0084-48	98-0084-78	98-0084-68	98-0084-88	98-0084-18
Large Adult Long (Burgundy)	31-40	320 x 158	810 x 183	98-0080-09	98-0084-39	98-0084-29	98-0080-39	98-0084-49	98-0084-79	98-0084-69	98-0084-89	98-0084-19
Thigh (Brown)	38-50	390 x 195	760 x 220	98-0080-10	98-0084-40	98-0084-30	98-0080-40	98-0084-50	98-0084-80	98-0084-70	98-0084-90	98-0084-20

Additional Specifications

Cuff Component	Material (No Latex, No PVC)
Cuff Shell	Nylon, 200 Denier
Hook	Nylon
Loop	Nylon
Thread	Nylon
Bladder	Nylon, 200 Denier
Bladder Tubing	TPE (thermoplastic elastomer)
Bladder Nipple	Polyurethane
Outer Label (If Applicable)	Polyester

Regulatory Requirements

Meets AHA, AMI-SP10 Guidelines

In compliance with Council Directive 93/42/EEC Medical Devices Directive

Meets ISO 10993 standard



Two-Tube Design Specifications



Cuff Size (Print Color)	Range (cm)	Bladder Length x Width (mm)	Overall Length x Width (mm)	Female Quick Connect	Male Quick Connect	Female & Male Quick Connects	Threaded (Screw)	Male Bayonet (Metal)	Male Bayonet (Plastic)	Female Slip Luer (Plastic)	No Connector
Infant (Orange)	8-13	110 x 38	230 x 63	98-0082-01	98-0085-31	98-0085-61	98-0085-21	98-0085-41	98-0085-71	98-0082-21	98-0085-11
Child (Green)	12-19	150 x 68	315 x 93	98-0082-02	98-0085-32	98-0085-62	98-0085-22	98-0085-42	98-0085-72	98-0082-22	98-0085-12
Child Long (Green)	12-19	150 x 68	445 x 93	98-0082-03	98-0085-33	98-0085-63	98-0085-23	98-0085-43	98-0085-73	98-0082-23	98-0085-13
Small Adult (Royal Blue)	17-25	193 x 93	440 x 118	98-0082-04	98-0085-34	98-0085-64	98-0085-24	98-0085-44	98-0085-74	98-0082-24	98-0085-14
Small Adult Long (Royal Blue)	17-25	193 x 93	515 x 118	98-0082-05	98-0085-35	98-0085-65	98-0085-25	98-0085-45	98-0085-75	98-0082-25	98-0085-15
Adult (Navy Blue)	23-33	265 x 125	505 x 150	98-0082-06	98-0085-36	98-0085-66	98-0085-26	98-0085-46	98-0085-76	98-0082-26	98-0085-16
Adult Long (Navy Blue)	23-33	265 x 125	645 x 150	98-0082-07	98-0085-37	98-0085-67	98-0085-27	98-0085-47	98-0085-77	98-0082-27	98-0085-17
Large Adult (Burgundy)	31-40	320 x 158	645 x 183	98-0082-08	98-0085-38	98-0085-68	98-0085-28	98-0085-48	98-0085-78	98-0082-28	98-0085-18
Large Adult Long (Burgundy)	31-40	320 x 158	810 x 183	98-0082-09	98-0085-39	98-0085-69	98-0085-29	98-0085-49	98-0085-79	98-0082-29	98-0085-19
Thigh (Brown)	38-50	390 x 195	760 x 220	98-0082-10	98-0085-40	98-0085-70	98-0085-30	98-0085-50	98-0085-80	98-0082-30	98-0085-20

Additional Specifications

Cuff Component	Material (No Latex, No PVC)
Cuff Shell	Nylon, 200 Denier
Hook	Nylon
Loop	Nylon
Thread	Nylon
Bladder	Nylon, 200 Denier
Bladder Tubing	TPE (thermoplastic elastomer)
Bladder Nipple	Polyurethane
Outer Label (If Applicable)	Polyester

Regulatory Requirements

Meets AHA, AMI-SP10 Guidelines

In compliance with Council Directive 93/42/EEC Medical Devices Directive

Meets ISO 10993 standard



SunTech Medical, Inc.
507 Airport Boulevard, Suite 117
Morrisville, NC 27560-8200 USA
Tel: + 1.919.654.2300
1.800.421.8626
Fax: + 1.919.654.2301

www.SunTechMed.com

SunTech Medical, Ltd.
Oakfield Estate
Eynsham, Oxfordshire
England
Tel: + 44 (0) 1865.884.234
Fax: + 44 (0) 1865.884.235

www.SunTechMed.com

SunTech Medical, Ltd.
Level 19, Two IFC
8, Finance Street, Central
Hong Kong
Tel: +852.2251.1949
Fax: +852.2251.1950

www.CN.SunTechMed.com

SunTech Medical, Ltd.
2/F of Building A, Jinxiongda Technology Park
Guanlan, Bao'an District Shenzhen, 518110, PRC
Tel: +86-755-2958 8810
+86-755-29588986 (Sales)
+86-755-29588665 (Service)
Fax: +86-755-2958 8829

www.CN.SunTechMed.com

part #82-0101-00 Rev. A