











| | | | | One Tube | | | | Two Tube | | | |
|------------------------------|------------|-----------------------------|-----------------------------|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|--------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|--------------|
| Design Specifications | | | |  |  |  | |  |  |  | |
| Cuff Size (Print Color) | Range (cm) | Bladder Length x Width (mm) | Overall Length x Width (mm) | Female Quick Connect | Threaded (Screw) | Male Bayonet (Plastic) | No Connector | Female Quick Connect | Threaded (Screw) | Male Bayonet (Plastic) | No Connector |
| Infant (Orange) | 8-13 | 120 x 54 | 228 x 59 | 98-0400-00 | 98-0400-20 | 98-0400-40 | 98-0400-60 | 98-0400-10 | 98-0400-30 | 98-0400-50 | 98-0400-70 |
| Child (Green) | 12-19 | 160 x 83 | 314 x 88 | 98-0400-01 | 98-0400-21 | 98-0400-41 | 98-0400-61 | 98-0400-11 | 98-0400-31 | 98-0400-51 | 98-0400-71 |
| Child Long (Green) | 12-19 | 160 x 83 | 440 x 88 | 98-0400-02 | 98-0400-22 | 98-0400-42 | 98-0400-62 | 98-0400-12 | 98-0400-32 | 98-0400-52 | 98-0400-72 |
| Small Adult (Turquoise) | 17-25 | 203 x 108 | 440 x 113 | 98-0400-03 | 98-0400-23 | 98-0400-43 | 98-0400-63 | 98-0400-13 | 98-0400-33 | 98-0400-53 | 98-0400-73 |
| Small Adult Long (Turquoise) | 17-25 | 203 x 108 | 515 x 113 | 98-0400-04 | 98-0400-24 | 98-0400-44 | 98-0400-64 | 98-0400-14 | 98-0400-34 | 98-0400-54 | 98-0400-74 |
| Adult (Navy Blue) | 23-33 | 275 x 140 | 505 x 145 | 98-0400-05 | 98-0400-25 | 98-0400-45 | 98-0400-65 | 98-0400-15 | 98-0400-35 | 98-0400-55 | 98-0400-75 |
| Adult Long (Navy Blue) | 23-33 | 275 x 140 | 640 x 145 | 98-0400-06 | 98-0400-26 | 98-0400-46 | 98-0400-66 | 98-0400-16 | 98-0400-36 | 98-0400-56 | 98-0400-76 |
| Large Adult (Burgundy) | 31-40 | 330 x 173 | 645 x 178 | 98-0400-07 | 98-0400-27 | 98-0400-47 | 98-0400-67 | 98-0400-17 | 98-0400-37 | 98-0400-57 | 98-0400-77 |
| Large Adult Long (Burgundy) | 31-40 | 330 x 173 | 815 x 178 | 98-0400-08 | 98-0400-28 | 98-0400-48 | 98-0400-68 | 98-0400-18 | 98-0400-38 | 98-0400-58 | 98-0400-78 |
| Thigh (Brown) | 38-50 | 400 x 210 | 755 x 215 | 98-0400-09 | 98-0400-29 | 98-0400-49 | 98-0400-69 | 98-0400-19 | 98-0400-39 | 98-0400-59 | 98-0400-79 |

| | | | | One Tube | | Two Tube | | | | | |
|-------------------------|------------|-----------------------------|-----------------------------|------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|--|--|--|--|
| Design Specifications | | | |  |  |  |  | | | | |
| Cuff Size (Print Color) | Range (cm) | Bladder Length x Width (mm) | Overall Length x Width (mm) | Male Slip Luer (Plastic) | Male Bayonet (Plastic) | Male Slip Luer (Plastic) | Male Bayonet (Plastic) | | | | |
| Neonate #1 (Orange) | 3-6 | 65 x 21 | 95 x 21 | 98-0400-80 | | 98-0400-85 | | | | | |
| Neonate #2 (Turquoise) | 4-8 | 75 x 28 | 113 x 28 | 98-0400-81 | | 98-0400-86 | | | | | |
| Neonate #3 (Green) | 6-11 | 95 x 36 | 140 x 36 | 98-0400-82 | | 98-0400-87 | | | | | |
| Neonate #4 (Navy Blue) | 7-13 | 121 x 44 | 175 x 44 | 98-0400-83 | | 98-0400-88 | | | | | |
| Neonate #5 (Burgundy) | 8-15 | 130 x 47 | 184 x 47 | 98-0400-84 | 98-0400-90 | 98-0400-89 | 98-0400-91 | | | | |

Adapters



Y Adapter, 2 Barb-to-barb (1/8") to 1 Metal Female Bayonet

98-0225-00



Y Adapter, 2 Barb-to-barb (1/8") to 1 Female Push-Button Bayonet

98-0226-00



Y Adapter, 2 Barb-to-barb (5/32") to 1 Metal Female Bayonet

98-0227-00



Y Adapter, 2 Barb-to-barb (5/32") to 1 Female Push-Button Bayonet

98-0228-00

Additional Specifications

| Cuff Component | Material |
|----------------|-----------------------------------------------------------------------------------|
| Cuff | Polyester (Dupont Sontara – Style 8001) with PVC (Polyvinyl chloride) - DEHP Free |
| Hook | PVC - DEHP Free |
| Loop | Nylon |
| Nipple | PVC - DEHP Free |
| Tubing | PVC - DEHP Free |



SunTech Medical, Inc.
507 Airport Boulevard, Suite 117
Morrisville, NC 27560-8200 USA
Tel: + 1.919.654.2300
1.800.421.8626
Fax: + 1.919.654.2301

www.SunTechMed.com

SunTech Medical, Ltd.
Oakfield Estate
Eynsham, Oxfordshire
England
Tel: + 44 (0) 1865.884.234
Fax: + 44 (0) 1865.884.235

www.SunTechMed.com

SunTech Medical, Ltd.
Level 19, Two IFC
8, Finance Street, Central
Hong Kong
Tel: +852.2251.1949
Fax: +852.2251.1950

www.CN.SunTechMed.com

SunTech Medical, Ltd.
2/F of Building A, Jinxiongda Technology Park
Guanlan, Bao'an District Shenzhen, 518110, PRC
Tel: +86-755-2958 8810
+86-755-29588986 (Sales)
+86-755-29588665 (Service)
Fax: +86-755-2958 8829

www.CN.SunTechMed.com

part #82-0110-00 Rev. A