



SunTech Vet20

Technical Specifications

The Vet20 Includes...

Vet20 System Monitor
 Veterinary BP Cuffs Sizes 1-6
 Patient Hose
 User guide
 Standard 2 year warranty

Vet20 System Monitors

99-0135-00 Vet20 Monitor



SunTech Veterinary Blood Pressure Cuffs

98-0400-80	Cuff #1	Box of 20	Range 3-8 cm
98-0400-81	Cuff #2	Box of 20	Range 4-8 cm
98-0400-82	Cuff #3	Box of 20	Range 6-11 cm
98-0400-83	Cuff #4	Box of 20	Range 7-13 cm
98-0400-84	Cuff #5	Box of 20	Range 8-15 cm
98-0400-F1	Cuff #6	Box of 20	Range 12-19 cm
98-0400-F3	Cuff #7	Box of 20	Range 17-25 cm
98-0240-00	Cuffs 1-6	Pack of 6	
91-0028-72	Patient Hose		

Accessories

39-0195-00 Vet20 Armour, Flamingo (Pink)
 39-0195-01 Vet20 Armour, Peacock (Blue)
 19-0018-00 AC Adapter

Extended Warranty Program

83-0021-00 One Year Extended Warranty
 83-0021-01 Three Years Extended Warranty

Extended Warranty Terms and Conditions

The Extended Warranty must be purchased either:

- At the time of the original sale
- Prior to the end of the initial warranty period
- Prior to the end of the purchased extended warranty period



General Specifications

SYS:	40 to 265 mmHg
DIA:	20 to 200 mmHg
MAP:	27 to 222 mmHg
HR:	25 to 300 mmHg
Device Weight:	1.3 lbs (with batteries) 590 grams (with batteries)
Dimensions:	6.25" x 5" x 5.25" 15.9 cm x 12.7 cm x 13.3 cm
Warranty:	2 years on monitor



507 Airport Boulevard, Suite 117
Morrisville, NC 27560-8200 USA
Tel: + 1.919.654.2300
1.800.421.8626
Fax: + 1.919.654.2301



Oakfield Industrial Estate
Stanton Harcourt Road
Eynsham, Oxfordshire OX29 4TS
England
Tel: + 44 (0) 1865.884.234
Fax: + 44 (0) 1865.884.235

105 HuanGuan South Road,
Suite 15 2~3/F
DaHe Community Guanlan,
LongHua District, Shenzhen
GuangDong PRC 518110
Tel.: + 86.755.29588810
+ 86.755.29588986 (Sales)
+ 86.755.29588665 (Service)
Fax: + 86.755.29588829
SunTechMed.com.CN

SunTechMed.com

SunTechMed.com


SunTech Medical[®]
The Difference in Clinical Grade[™]

