



 SunTech[®] **Vet25[™]**
BP MONITOR



Spot-Check and Continuous BP Monitor
for Companion Animals

SunTech Vet25

The portable and rechargeable SunTech® Vet25™ monitor quickly and easily provides reliable blood pressure (BP) measurements for wellness exams, spot-checks, and while monitoring companion animals pre- and post-procedure. Intervals can be set to automatically take BP readings during a monitoring period. Utilizing SunTech's trusted Advantage™ VET BP technology, the Vet25 quietly performs motion-tolerant BP measurements, delivering accurate results. The intuitive, touchscreen interface provides users all of the necessary features for reviewing and analyzing BP data.

Standard Features

Reliable:

Due to SunTech's motion tolerant technology, the Vet25 has a high success rate for BP readings.

Simple to Use:

Intuitive touch screen interface provides all necessary features for reviewing, analyzing, and transferring BP data. Designed for animal hospitals, not researchers.

Accurate:

BP algorithm was specifically developed for companion animals by experts at a leading vet school.

Continuous Monitoring: Monitoring during procedures is simple with user-defined intervals for BP measurement.

Memory:

Up to 960 data points are retained, even when device is powered off. Users can select specific measurements from device memory and perform automated BP averaging.

Portable:


Battery-operated device easily moved to different monitoring locations.

Rechargeable:

Battery recharges when device is plugged into AC power, reducing disposable battery waste.

Data Transfer:

Data can be transferred to a data trending and reporting software application for a PC using Bluetooth. Use the PC application to export reports in either PDF or CSV format.





Technical Specifications

The Vet25 Includes:

- Vet25 System Monitor
- Protective Armour in Flamingo Pink, Peacock Blue, or Tree Frog Green
- Veterinary Blood Pressure Cuffs Sizes 1-6
- Patient Hose, 6 ft (1.8 m)
- AC Adapter

General Specifications

Systolic:	40 to 265 mmHg
Diastolic:	20 to 200 mmHg
Mean Arterial Pressure:	27 to 222 mmHg
Heart Rate:	25 to 300 BPM
Device Weight:	2.1 lbs (with battery), 953 grams (with battery)
Dimensions:	6.25" x 5" x 5.25", 15.9 cm x 12.7 cm x 13.3 cm
Warranty:	2 years on monitor 90 days on cuffs/hose

Accessories

SunTech Vet25 Armour

The removable protective armour is available in 3 colors and features a convenient handle.



Flamingo Pink



Peacock Blue



Tree Frog Green

AC Adapter

The Vet25 is a portable device that can easily travel from room to room when needed. Simply plug in the AC Adapter to recharge the lithium-ion battery in the Vet25. The AC Adapter includes 4 plug configurations to meet most international power requirements.

SunTech Veterinary Cuffs



Soft materials and rounded corners provide a gentle cuff solution. SunTech Medical BP cuffs are color-coded for quick size selection and range indicators ensure accurate cuff sizing and fit. Included with the Vet25 BP device are cuffs size 1–6. **SunTech Vet BP cuffs can be used with other veterinary BP monitors! Ask about our SunTech Vet Cuff Conversion Kit.**


SunTech Medical®
The Difference in Clinical Grade™

With more than 30 years dedicated to developing our world-class blood pressure technology, we are confident that we don't miss a beat - so you won't either. It's our legacy. It's SunTech Medical. **It's The Difference In Clinical Grade.**



SunTech Medical, Inc.
507 Airport Boulevard, Suite 117
Morrisville, NC 27560-8200 USA
Tel: + 1.919.654.2300
1.800.421.8626
Fax: +1.919.654.2301
SunTechMed.com

EC REP

SunTech Medical, Ltd.
Oakfield Industrial Estate
Stanton Harcourt Road
Eynsham, Oxfordshire OX29 4TS
England
Tel: + 44 (0) 1865.884.234
Fax: + 44 (0) 1865.884.235
SunTechMed.com

SunTech Medical (Shenzhen) Co., Ltd.
105 HuanGuan South Road, Suite 15 2~3/F
DaHe Community Guanlan,
LongHua District, Shenzhen
GuangDong PRC 518110
Tel.: + 86.755.29588810
+ 86.755.29588986 (Sales)
+ 86.755.29588665 (Service)
Fax: + 86.755.29588829
SunTechMed.com.CN





SunTech Vet25 Part Numbers

Vet25 System Monitors

99-0170-00	Vet25 Monitor, Flamingo Pink Armour
99-0170-01	Vet25 Monitor, Peacock Blue Armour
99-0170-02	Vet25 Monitor, Tree Frog Green Armour

SunTech Veterinary BP Cuffs

98-0400-80	Cuff #1	Box of 20	Range 3-6 cm
98-0400-81	Cuff #2	Box of 20	Range 4-8 cm
98-0400-82	Cuff #3	Box of 20	Range 6-11 cm
98-0400-83	Cuff #4	Box of 20	Range 7-13 cm
98-0400-84	Cuff #5	Box of 20	Range 8-15 cm
98-0400-F1	Cuff #6	Box of 20	Range 12-19 cm
98-0400-F3	Cuff #7	Box of 20	Range 17-25 cm
98-0240-00	Cuffs 1-6	Pack of 6	
91-0028-75	Patient Hose, 6 ft (1.8 m)		
91-0028-76	Patient Hose, 10 ft (3 m)		

Accessories

39-0195-00	Flamingo Pink Armour
39-0195-01	Peacock Blue Armour
39-0195-02	Tree Frog Green Armour
19-0020-00	AC Adapter
98-0900-01	Lithium-ion Battery Replacement

Extended Warranty Program

83-0022-00	One Year Extended Warranty
83-0022-01	Three Years Extended Warranty

Extended Warranty Terms and Conditions

The Extended Warranty can be purchased either:

- At the time of the original sale
- Prior to the end of the initial warranty period