



 SunTech® **Vet25™**  
BP MONITOR



**Spot-Check and Continuous BP Monitor**  
for Companion Animals

# SunTech Vet25

The portable and rechargeable SunTech® Vet25™ monitor quickly and easily provides reliable blood pressure (BP) measurements for wellness exams, spot-checks, and while monitoring companion animals pre- and post-procedure. Intervals can be set to automatically take BP readings during a monitoring period. Utilizing SunTech's trusted Advantage™ VET BP technology, the Vet25 quietly performs motion-tolerant BP measurements, delivering accurate results. The intuitive, touchscreen interface provides users all of the necessary features for reviewing and analyzing BP data.

## Standard Features

**Reliable:** Due to SunTech's motion tolerant technology, the Vet25 has a high success rate for BP readings.

**Simple to Use:** Intuitive touch screen interface provides all necessary features for reviewing, analyzing, and transferring BP data. Designed for animal hospitals, not researchers.

**Accurate:** BP algorithm was specifically developed for companion animals by experts at a leading vet school.


**Continuous Monitoring:** Monitoring during procedures is simple with user-defined intervals for BP measurement.

**Memory:** Up to 960 data points are retained, even when device is powered off. Users can select specific measurements from the device memory and perform automated BP averaging.

**Portable:** Battery-operated device easily moved to different monitoring locations.

**Rechargeable:** Battery recharges when device is plugged into AC power, reducing disposable battery waste.

**Data Transfer:** Data can be transferred to a data trending and reporting software application for a PC using Bluetooth. Use the PC application to export reports in either PDF or CSV format.





## Technical Specifications

### The Vet25 Includes:

- Vet25 System Monitor
- Protective Armour in Flamingo Pink, Peacock Blue, or Tree Frog Green
- Veterinary Blood Pressure Cuffs Sizes 1-6
- Patient Hose, 6 ft (1.8 m)
- AC Adapter

### General Specifications

<b>Systolic:</b>	40 to 265 mmHg
<b>Diastolic:</b>	20 to 200 mmHg
<b>Mean Arterial Pressure:</b>	27 to 222 mmHg
<b>Heart Rate:</b>	25 to 300 BPM
<b>Device Weight:</b>	2.1 lbs (with battery), 953 grams (with battery)
<b>Dimensions:</b>	6.25" x 5" x 5.25", 15.9 cm x 12.7 cm x 13.3 cm
<b>Warranty:</b>	2 years on monitor 90 days on cuffs/hose

## Accessories

### SunTech Vet25 Armour

The removable protective armour is available in 3 colors and features a convenient handle.



*Flamingo Pink*



*Peacock Blue*



*Tree Frog Green*

### AC Adapter

The Vet25 is a portable device that can easily travel from room to room when needed. Simply plug in the AC Adapter to recharge the lithium-ion battery in the Vet25. The AC Adapter includes 4 plug configurations to meet most international power requirements.

## SunTech Veterinary Cuffs



Soft materials and rounded corners provide a gentle cuff solution. SunTech Medical BP cuffs are color-coded for quick size selection and range indicators ensure accurate cuff sizing and fit. Included with the Vet25 BP device are cuffs size 1–6.

**SunTech Vet BP cuffs can be used with other veterinary BP monitors! Ask about our SunTech Vet Cuff Conversion Kit.**

  
**SunTech Medical**<sup>®</sup>  
The Difference in Clinical Grade<sup>™</sup>

With more than 30 years dedicated to developing our world-class blood pressure technology, we are confident that we don't miss a beat - so you won't either. It's our legacy. It's SunTech Medical. ***It's The Difference In Clinical Grade.***



SunTech Medical, Inc.  
507 Airport Boulevard, Suite 117  
Morrisville, NC 27560-8200 USA  
Tel: + 1.919.654.2300  
1.800.421.8626  
Fax: +1.919.654.2301  
SunTechMed.com

EC REP

SunTech Medical, Ltd.  
Oakfield Industrial Estate  
Stanton Harcourt Road  
Eynsham, Oxfordshire OX29 4TS  
England  
Tel: + 44 (0) 1865.884.234  
Fax: + 44 (0) 1865.884.235  
SunTechMed.com

SunTech Medical (Shenzhen) Co., Ltd.  
105 HuanGuan South Road, Suite 15 2~3/F  
DaHe Community Guanlan,  
LongHua District, Shenzhen  
GuangDong PRC 518110  
Tel.: + 86.755.29588810  
+ 86.755.29588986 (Sales)  
+ 86.755.29588665 (Service)  
Fax: + 86.755.29588829  
SunTechMed.com.CN





# SunTech Vet25 Part Numbers

## Vet25 System Monitors

- 99-0170-00 Vet25 Monitor, Flamingo Pink Armour
- 99-0170-01 Vet25 Monitor, Peacock Blue Armour
- 99-0170-02 Vet25 Monitor, Tree Frog Green Armour

## SunTech Veterinary BP Cuffs

- |            |                            |           |                |
|------------|----------------------------|-----------|----------------|
| 98-0400-80 | Cuff #1                    | Box of 20 | Range 3-6 cm   |
| 98-0400-81 | Cuff #2                    | Box of 20 | Range 4-8 cm   |
| 98-0400-82 | Cuff #3                    | Box of 20 | Range 6-11 cm  |
| 98-0400-83 | Cuff #4                    | Box of 20 | Range 7-13 cm  |
| 98-0400-84 | Cuff #5                    | Box of 20 | Range 8-15 cm  |
| 98-0400-F1 | Cuff #6                    | Box of 20 | Range 12-19 cm |
| 98-0400-F3 | Cuff #7                    | Box of 20 | Range 17-25 cm |
| 98-0240-00 | Cuffs 1-6                  | Pack of 6 |                |
| 91-0028-75 | Patient Hose, 6 ft (1.8 m) |           |                |
| 91-0028-76 | Patient Hose, 10 ft (3 m)  |           |                |

## Accessories

- 39-0195-03 Flamingo Pink Armour
- 39-0195-04 Peacock Blue Armour
- 39-0195-02 Tree Frog Green Armour
- 19-0020-00 AC Adapter
- 98-0900-01 Lithium-ion Battery Replacement

## Extended Warranty Program

- 83-0022-00 One Year Extended Warranty
- 83-0022-01 Three Years Extended Warranty

## Extended Warranty Terms and Conditions

The Extended Warranty can be purchased either:

- At the time of the original sale
- Prior to the end of the initial warranty period